

# norma

Diagnostic Solutions

## icon 3

### A Smart Hematology Analyzer

- Smallest sample volume makes it ideal for pediatric use
- Smallest reagent consumption
- A4 paper size footprint
- ER mode for use at the point-of-care
- NormaAccess
  - Remote support tool



[www.normadiagnostika.com](http://www.normadiagnostika.com)

## TECHNICAL DATA

Parameters	WBC, LYM, MID, GRA, LYM%, MID%, GRA%, RBC, HGB, MCV, HCT, MCH, MCHC, RDWcv, RDWsd, PLT, MPV, PCT, PDWcw, PDWsd, P-LCR%, P-LCC
Measurement technology	Volumetrical impedance method for WBC and RBC Photometrical method for HGB
Histograms	WBC, RBC, PLT
Sampling mode	Closed and open vial mode
Aspirated sample volume	Closed vial: 38 µl Open vial: 13 µl
Processed sample volume	2.5 µl
Throughput	60 tests per hour
Precision (CV%)	WBC < 3%, RBC < 2%, PLT < 5% MCV < 1%, HGB < 2%
Storage capacity	40,000 test results incl. histograms
User interface	LCD, 10.1", 1280 × 800 with capacitive touchscreen, portrait
Connectivity	2 x USB ports, Ethernet, LIS (HL7 and ASTM)
Dimensions (HxWxD)	280 x 216 x 320 mm
Net weight	8.7 kg

Norma Icon is ideal for doctors' offices, emergency departments, small to medium size labs and polyclinics.



- Microfluidic technology
- Patented, innovative solutions
- Trendsetting design



Shortlisted for the Senetics Innovation Award 2014  
Category "Most Innovative Product in the Health Care Sector"

reddot award 2014  
winner

icon

Norma Instruments Zrt.

Papírgyár u. 58-59., 1038 Budapest  
Hungary

Phone: +36 1 815 4370

E-mail: office@normadiagnostika.com



SOWAR PRIVATE LIMITED

CSC, C-9, VASANT KUNJ, NEW DELHI - 110070

Tel: 011 - 41767519 - 22, Fax: 011 - 26122809

Email: sales@sowarmedical.com; sowar@sowarmedical.com

www.sowarmedical.com

### Sales Offices:

Chennai, Bangalore, Hyderabad, Kerala  
Mumbai, Pune, Nagpur, Ahmedabad  
Kolkata, Orissa, Guwahati, Bihar  
Chhattisgarh, Madhya Pradesh  
Rajasthan, Punjab, Banaras